Dr.Salmon Newsletter

# World News & Views

## -Letters from Dr.Salmon, NSU-

Dear readers,

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**Thomas O. Salmon, OD, PhD, FAAO**  
Professor, Northeastern State University

**VIA AIR MAIL**

CooperVision **4e** Program  
enhance each and every contact lens experience.

## ARVO Reviews

VISIONARY GENOMICS  
FORT LAUDERDALE, FL  
MAY 1-5, 2011

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Influence of Commercial Soft Multifocal Contact Lenses on Peripheral Refraction and Aberrations

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Evidence of a Compensation Mechanism between Corneal and Internal Ocular Optical Aberrations in Patients with Keratoconus

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**On-Road Driving Performance in Moderate and Advanced Glaucoma Patients**

**演者:** Victoria H. Yom (Washington University School of Medicine, St. Louis, Missouri, USA)

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**Yoga As A Treatment Modality For Persons With Severe Visual Impairment: A Feasibility Study**

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**Corneal Topographic Changes During Continuous Silicone Hydrogel Contact Lens Wear**

**演者:** Scott W. Yeates (Case Western Reserve University and University Hospitals Eye Institute, Cleveland, Ohio, USA)

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**Evaluating the Biocidal Efficacy of Multi-purpose and Hydrogen Peroxide Solutions Against Fungal Isolates Associated with Ocular Infection**

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**The Effect of Lens Modulus on Insertion Comfort with Silicone Hydrogel Lenses**

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
演者: Junko Tsukiyama (Kinki University School of Medicine, Osaka-Sayama City, Japan)

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**5533 Contamination of Silicone Hydrogel Contact Lenses by Cosmetics and the Cleaning Effect of Contact Lens Daily Cleaners**  
 Junko Tsukiyama<sup>1</sup>, Yoko Miyamoto<sup>1</sup>, Masahiko Fukuda<sup>1</sup>, Yoshikazu Shimomura<sup>1</sup>, Hirohiko Minori<sup>1</sup>, Jiro Tsukiyagi<sup>2</sup>  
 Department of Ophthalmology, Kinki University Faculty of Medicine, Osaka-Sayama city, Japan<sup>1</sup>, CIBA Vision K.K., Tokyo, Japan<sup>2</sup>

**Purpose:** We investigated the cleaning effect of contact lens daily cleaners on silicone hydrogel contact lenses contaminated by cosmetics. We used two types of daily cleaners: MPS (MPS, Ciba Vision, Tokyo, Japan) and Miraflow (Miraflow, Ciba Vision, Tokyo, Japan). We used two types of contact lenses: senofilcon A (senofilcon A, Ciba Vision, Tokyo, Japan) and balafilcon A (balafilcon A, Ciba Vision, Tokyo, Japan). We used two types of cosmetics: pencil eyeliner (pencil eyeliner, Ciba Vision, Tokyo, Japan) and cleansing oil (cleansing oil, Ciba Vision, Tokyo, Japan).

**Materials:** senofilcon A, balafilcon A, asmofilcon A, pencil eyeliner, cleansing oil, MPS, Miraflow.

**Methods:** We used two types of contact lenses: senofilcon A and balafilcon A. We used two types of cosmetics: pencil eyeliner and cleansing oil. We used two types of daily cleaners: MPS and Miraflow. We used two types of solutions: physiological saline and distilled water.

**Results:** We investigated the cleaning effect of contact lens daily cleaners on silicone hydrogel contact lenses contaminated by cosmetics. We used two types of daily cleaners: MPS and Miraflow. We used two types of contact lenses: senofilcon A and balafilcon A. We used two types of cosmetics: pencil eyeliner and cleansing oil. We used two types of solutions: physiological saline and distilled water.

**Conclusion:** Using a 30-second saline wash and 20-second rubbing with Miraflow is most effective in removing cosmetic contamination from lenses, but it could not remove all contamination. Contact lens wearers who use eye cosmetics should be careful to avoid contaminating their lenses with cosmetics.

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Comparison of Morning and Afternoon Osmolarity in Silicone Hydrogel Wearers and Non-contact Lens Wearers

演者: Ryan Dimit (University of Houston, Houston, TX)

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Optical Modeling And Analysis Of Peripheral Optics Of Contact Lenses

演者: Jie Shen (Western University of Health Sciences, Pomona, CA)

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**Optical Modeling And Analysis Of Peripheral Optics Of Contact Lenses** Program 6547

Jie Shen<sup>1</sup>, Frank Spora<sup>2</sup>, Larry Thibos<sup>2</sup>  
<sup>1</sup>Western University of Health Sciences, College of Optometry, Pomona, CA  
<sup>2</sup>School of Optometry, Indiana University, Bloomington, IN

**Abstract**

Background: Contact lenses (CL) are widely used to correct refractive errors. The purpose of this study was to use optical modeling to evaluate peripheral optics of various contact lenses in various. Zemax EE optical design software was used to model peripheral optics of rigid gas permeable lenses (RGPCL) and soft contact lenses (SCL). Profiles of lower and higher order aberrations were compared and 40° eccentricity respectively on 10° eccentricity. Results: The results have investigated peripheral optics of CL in relation to the study. We used optical design software to evaluate the peripheral optics of left and right CL in relation to better understand their complex optical impact on central image formation.

**Introduction**

Types with peripheral hyperopic refractive error or degraded peripheral vision are at risk for developing myopia. Contact lenses are widely used for correction of central visual refractive errors, thereby diffuse and uniform. Their effect on peripheral vision has been studied in both experimental and clinical settings. Results: The results have investigated peripheral optics of CL in relation to the study. We used optical design software to evaluate the peripheral optics of left and right CL in relation to better understand their complex optical impact on central image formation.

**Methods**

Zemax EE optical design software was used to model peripheral optics of RGPCL (TET design, Menicon) and SCL (Acuvue 2, Johnson & Johnson). Aberrations profiles (up to 10° eccentricity) were compared on 40° steps up to 40° off-axis eccentricity along the horizontal meridian for CL with dioptric power from -3.0 D to 0 D. Glare (ringing) results were analyzed and compared with published experimental data.<sup>1</sup> Corneal aberrations data were simulated using the following equation:<sup>2</sup>

$$W_{040} = \frac{1}{2} C_4 r^4 + \frac{1}{2} C_6 r^6 + \frac{1}{2} C_8 r^8$$

**Results**

Fig. 1: Both CL types produced less negative dioptric power profiles. As off-axis eccentricity increased for RGPCL, so did SCL, but SCL was 2.0 D.

Fig. 2: CL induced "against-the-RSA" astigmatism in the periphery along the horizontal meridian. RGPCL, manufacturers, J.

Fig. 3: A comparison of optical modeling results with experimental data showed good agreement within a field angle of 20°.

Fig. 4: Both CL types showed increasing positive value of peripheral relative spherical aberration (PRSA) when the central dioptric power varied from -3.0 D to 0.0 D. SCL also increased systematically in the periphery with increasing central power at any of the CL types.

Fig. 5: Optical modeling of major higher-order aberrations (spherical, RGPCL and SCL) with smaller central power revealed the most increased amount of coma and negative spherical aberration toward the periphery. The amount of these higher order aberrations increased with increased central focal power for both lens types.

**Conclusion**

Compared to SCL, RGPCL with equal central focal power produced less spherical aberration in the periphery. Both lens types produced comparable amounts of astigmatism and major higher-order aberrations with increased field angle, in particular coma and spherical aberration. All lens induced effects depended on the central dioptric power and increased when the latter increased. Optical modeling results of CL revealed approximately with previously obtained experimental findings on the change of PRSA with CL correction. Thus the modeling of peripheral optics of CL in relation can be used to predict field-of-vision changes after CL correction in myopic eyes.

**References**

1. Sheipholi, J., Basset, F., et al. (2012). "Against-the-RSA" astigmatism in young patients. *Ophthalmologica* 125(1): 299-313.
2. Williams, C. J. (1987). "Using wavefront measurements for the detection and classification of clinical practice." *Ophthalmic Physiol Opt* 7(5): 274-80.
3. Mann, D. G., Sheiko, R. J., et al. (2009). "Peripheral refractive and ocular shape in children." *Invest Ophthalmol Vis Sci* 50(12): 3222-30.
4. Schachar, R. S. (2004). "Astigmatism response to myopia in children." *Optom* 95(10): 22-25.
5. Williams, J. D., et al. (2010). "Measurement of Eye Growth and the Role of Myopia." *Optom* 91(1): 16-24.
6. Sheiko, R., et al. (2010). "Peripheral refractive aberrations produce ocular elongation in chicks." *Invest Ophthalmol Vis Sci* 51(12): 6149.
7. Dainoff, A., Schaeffel, F. (2007). "Local changes in eye growth induced by peripheral focal refractive error despite central accommodation." *Visum Res* 47(8): 879-88.
8. Williams, M., Schaeffel, F. (2006). "A simple mechanism for accommodation without eye size accommodation in chicks." *Visum Res* 46(3): 874-8.
9. Ponce, L., Sheiko, R., et al. (2010). "The effect of central lens power profile on peripheral refractive error." *AAO meeting* 2010. San Francisco. Program 210108.
10. Shen, J., Chak, A., et al. (2010). "Peripheral Refraction With and Without Contact Lens Correction." *Optom Vis Sci* 87(1): 60-65.
11. Thibos, L. N., Whitley, W., Horner, D. "Power spectra: an application of Fourier analysis to the description and statistical analysis of refractive error." *Optom Vis Sci* 10(7): 567-75.

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