



CooperVision®

頂点間距離補正表

眼鏡からコンタクトレンズの
パワーを求める換算表

1.5
2.0
2.5
3.0
3.5
4.0
4.5
5.0
5.5
6.0
6.5
7.0
7.5
8.0
11.0
12.0
13.0

| 眼鏡レンズ | 眼前 12mm | | 眼前 14mm | |
|-------|---------|-------|---------|-------|
| | - | + | - | + |
| 3.00 | 3.00 | 3.00 | 3.00 | 3.25 |
| .25 | 3.25 | 3.50 | 3.00 | 3.50 |
| .50 | 3.25 | 3.75 | 3.25 | 3.75 |
| .75 | 3.50 | 4.00 | 3.50 | 4.00 |
| 4.00 | 3.75 | 4.25 | 3.75 | 4.25 |
| .25 | 4.00 | 4.50 | 4.00 | 4.50 |
| .50 | 4.25 | 4.75 | 4.25 | 4.75 |
| .75 | 4.50 | 5.00 | 4.50 | 5.00 |
| 5.00 | 4.75 | 5.25 | 4.75 | 5.50 |
| .25 | 5.00 | 5.50 | 5.00 | 5.75 |
| .50 | 5.25 | 6.00 | 5.00 | 6.00 |
| .75 | 5.50 | 6.25 | 5.25 | 6.25 |
| 6.00 | 5.50 | 6.50 | 5.50 | 6.50 |
| .25 | 5.75 | 6.75 | 5.75 | 6.75 |
| .50 | 6.00 | 7.00 | 6.00 | 7.25 |
| .75 | 6.25 | 7.25 | 6.25 | 7.50 |
| 7.00 | 6.50 | 7.75 | 6.50 | 7.75 |
| .25 | 6.75 | 8.00 | 6.50 | 8.00 |
| .50 | 7.00 | 8.25 | 6.75 | 8.50 |
| .75 | 7.00 | 8.50 | 7.00 | 8.75 |
| 8.00 | 7.25 | 8.75 | 7.25 | 9.00 |
| .25 | 7.50 | 9.25 | 7.50 | 9.25 |
| .50 | 7.75 | 9.50 | 7.50 | 9.75 |
| .75 | 8.00 | 9.75 | 7.75 | 10.00 |
| 9.00 | 8.00 | 10.00 | 8.00 | 10.25 |
| .25 | 8.25 | 10.50 | 8.25 | 10.75 |
| .50 | 8.50 | 10.75 | 8.50 | 11.00 |
| .75 | 8.75 | 11.00 | 8.50 | 11.25 |
| 10.00 | 9.00 | 11.25 | 8.75 | 11.75 |
| .50 | 9.25 | 12.00 | 9.25 | 12.25 |
| 11.00 | 9.75 | 12.75 | 9.50 | 13.00 |
| .50 | 10.00 | 13.25 | 10.00 | 13.75 |
| 12.00 | 10.50 | 14.00 | 10.25 | 14.50 |
| .50 | 10.75 | 14.75 | 10.75 | 15.25 |
| 13.00 | 11.25 | 15.50 | 11.00 | 16.00 |
| .50 | 11.50 | 16.00 | 11.25 | 16.75 |
| 14.00 | 12.00 | 16.75 | 11.75 | 17.50 |
| .50 | 12.25 | 17.50 | 12.00 | 18.25 |
| 15.00 | 12.75 | 18.25 | 12.50 | 19.00 |
| .50 | 13.00 | 19.00 | 12.75 | 19.75 |
| 16.00 | 13.50 | 19.75 | 13.00 | 20.50 |
| .50 | 13.75 | 20.50 | 13.50 | 21.50 |
| 17.00 | 14.00 | 21.25 | 13.75 | 22.25 |
| .50 | 14.50 | 22.25 | 14.00 | 23.25 |
| 18.00 | 14.75 | 23.00 | 14.50 | 24.00 |
| .50 | 15.25 | 23.75 | 14.75 | 25.00 |
| 19.00 | 15.50 | 24.50 | 15.00 | 26.00 |
| .50 | 15.75 | 25.50 | 15.25 | 26.75 |
| 20.00 | 16.25 | 26.25 | 15.75 | 27.75 |

0
10
20
30
40
50
60
70
80
90
100
110
120
130
140
150