


度数選択早見表

眼鏡矯正球面度数	Cyl.-0.75		Cyl.-1.00		Cyl.-1.25		Cyl.-1.50		Cyl.-1.75		Cyl.-2.00		Cyl.-2.25		Cyl.-2.50		Cyl.-2.75		Cyl.-3.00	
	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数	球面度数	円柱度数
+5.00	+5.50		+5.25		+5.25		+5.25		+5.25		+5.25		+5.25		+5.25		+5.00		+5.00	
+4.75	+5.00		+5.00		+5.00		+5.00		+5.00		+5.00		+5.00		+5.00		+4.75		+4.75	
+4.50	+4.75		+4.75		+4.75		+4.75		+4.75		+4.75		+4.75		+4.75		+4.50		+4.50	
+4.25	+4.50		+4.50		+4.50		+4.50		+4.50		+4.50		+4.50		+4.50		+4.25		+4.25	
+4.00	+4.25		+4.25		+4.25		+4.25		+4.25		+4.25		+4.25		+4.00		+4.00		+4.00	
+3.75	+4.00		+4.00		+4.00		+3.75		+4.00		+3.75		+4.00		+3.75		+3.75		+3.75	
+3.50	+3.75		+3.50		+3.75		+3.50		+3.75		+3.50		+3.75		+3.50		+3.50		+3.50	
+3.25	+3.50		+3.25		+3.25		+3.25		+3.50		+3.25		+3.50		+3.25		+3.25		+3.25	
+3.00	+3.25		+3.00		+3.25		+3.00		+3.25		+3.00		+3.25		+3.00		+3.00		+3.00	
+2.75	+3.00		+2.75		+3.00	-1.25	+2.75		+3.00	-1.75	+2.75		+3.00	-2.25	+2.75		+2.75		+2.75	
+2.50	+2.75		+2.50		+2.75		+2.50		+2.75		+2.50		+2.75		+2.50		+2.50		+2.50	
+2.25	+2.50		+2.25		+2.50		+2.25		+2.50		+2.25		+2.50		+2.25		+2.25		+2.25	
+2.00	+2.25		+2.00		+2.25		+2.00		+2.25		+2.00		+2.25		+2.00		+2.00		+2.00	
+1.75	+2.00		+1.75		+2.00		+1.75		+2.00		+1.75		+2.00		+1.75		+1.75		+1.75	
+1.50	+1.75		+1.50		+1.75		+1.50		+1.75		+1.50		+1.75		+1.50		+1.50		+1.50	
+1.25	+1.50		+1.25		+1.50		+1.25		+1.50		+1.25		+1.50		+1.25		+1.25		+1.25	
+1.00	+1.25		+1.00		+1.25		+1.00		+1.25		+1.00		+1.25		+1.00		+1.00		+1.00	
+0.75	+0.75		+0.50		+0.75		+0.75		+0.75		+0.75		+0.75		+0.75		+0.75		+0.75	
+0.50	+0.50		+0.50		+0.50		+0.50		+0.50		+0.50		+0.50		+0.50		+0.50		+0.50	
+0.25	+0.25		+0.25		+0.25		+0.25		+0.25		+0.25		+0.25		+0.25		+0.25		+0.25	
±0.00	±0.00		±0.00		±0.00		±0.00		±0.00		±0.00		±0.00		±0.00		±0.00		±0.00	
-0.25	-0.25		-0.25		-0.25		-0.25		-0.25		-0.25		-0.25		-0.25		-0.25		-0.25	
-0.50	-0.50		-0.50		-0.50		-0.50		-0.50		-0.50		-0.50		-0.50		-0.50		-0.50	
-0.75	-0.75		-0.75		-0.75		-0.75		-0.75		-0.75		-0.75		-0.75		-0.75		-0.75	
-1.00	-1.00		-1.00		-1.00		-1.00		-1.00		-1.00		-1.00		-1.00		-1.00		-1.00	
-1.25	-1.25		-1.25		-1.25		-1.25		-1.25		-1.25		-1.25		-1.25		-1.25		-1.25	
-1.50	-1.50		-1.50		-1.50		-1.50		-1.50		-1.50		-1.50		-1.50		-1.50		-1.50	
-1.75	-1.75		-1.75		-1.75		-1.75		-1.75		-1.75		-1.75		-1.75		-1.75		-1.75	
-2.00	-1.75	-0.75	-2.00	-0.75	-2.00	-0.75	-2.00	-1.25	-2.00	-1.25	-2.00	-1.75	-2.00	-1.75	-2.00	-2.25	-2.00	-2.25	-2.00	-2.25
-2.25	-2.00		-2.25		-2.25		-2.25		-2.25		-2.25		-2.25		-2.25		-2.25		-2.25	
-2.50	-2.25		-2.50		-2.50		-2.50		-2.50		-2.50		-2.50		-2.50		-2.50		-2.50	
-2.75	-2.50		-2.75		-2.75		-2.75		-2.75		-2.75		-2.75		-2.75		-2.75		-2.75	
-3.00	-2.75		-3.00		-3.00		-3.00		-3.00		-2.75		-3.00		-2.75		-3.00		-3.00	
-3.25	-3.00		-3.25		-3.25		-3.00		-3.25		-3.00		-3.25		-3.00		-3.25		-3.25	
-3.50	-3.25		-3.50		-3.50		-3.25		-3.50		-3.25		-3.50		-3.25		-3.50		-3.50	
-3.75	-3.50		-3.75		-3.75		-3.50		-3.75		-3.50		-3.75		-3.50		-3.75		-3.75	
-4.00	-3.75		-4.00		-4.00		-3.75		-4.00		-3.75		-4.00		-3.75		-4.00		-4.00	
-4.25	-4.00		-4.25		-4.25		-4.00		-4.25		-4.00		-4.25		-4.00		-4.25		-4.25	
-4.50	-4.25		-4.50		-4.50		-4.25		-4.50		-4.25		-4.50		-4.25		-4.50		-4.50	
-4.75	-4.50		-4.75		-4.75		-4.50		-4.75		-4.50		-4.75		-4.50		-4.75		-4.75	
-5.00	-4.75		-5.00		-5.00		-4.75		-5.00		-4.75		-5.00		-4.75		-5.00		-5.00	
-5.25	-5.00		-5.25		-5.25		-5.00		-5.25		-5.00		-5.25		-5.00		-5.25		-5.25	
-5.50	-5.25		-5.50		-5.50		-5.25		-5.50		-5.25		-5.50		-5.25		-5.50		-5.50	
-5.75	-5.50		-5.75		-5.75		-5.50		-5.75		-5.50		-5.75		-5.50		-5.75		-5.75	
-6.00	-5.75		-6.00		-6.00		-5.75		-6.00		-5.75		-6.00		-5.75		-6.00		-6.00	
-6.25	-6.00		-6.25		-6.25		-6.00		-6.25		-6.00		-6.25		-6.00		-6.25		-6.25	
-6.50	-6.25		-6.50		-6.50		-6.25		-6.50		-6.25		-6.50		-6.25		-6.50		-6.50	
-6.75	-6.50		-6.75		-6.75		-6.50		-6.75		-6.50		-6.75		-6.50		-6.75		-6.75	
-7.00	-6.75		-7.00		-7.00		-6.75		-7.00		-6.75		-7.00		-6.75		-7.00		-7.00	
-7.25	-7.00		-7.25		-7.25		-7.00		-7.25		-7.00		-7.25		-7.00		-7.25		-7.25	
-7.50	-7.25		-7.50		-7.50		-7.25		-7.50		-7.25		-7.50		-7.25		-7.50		-7.50	
-7.75	-7.50		-7.75		-7.75		-7.50		-7.75		-7.50		-7.75		-7.50		-7.75		-7.75	
-8.00	-7.75		-8.00		-8.00		-7.75		-8.00		-7.75		-8.00		-7.75		-8.00		-8.00	
-8.25	-8.00		-8.25		-8.25		-8.00		-8.25		-8.00		-8.25		-8.00		-8.25		-8.25	
-8.50	-8.25		-8.50		-8.50		-8.25		-8.50		-8.25		-8.50		-8.25		-8.50		-8.50	
-8.75	-8.50		-8.75		-8.75		-8.50		-8.75		-8.50		-8.75		-8.50		-8.75		-8.75	
-9.00	-8.75		-9.00		-9.00		-8.75		-9.00		-8.75		-9.00		-8.75		-9.00		-9.00	
-9.25	-9.00		-9.25		-9.25		-9.00		-9.25		-9.00		-9.25		-9.00		-9.25		-9.25	
-9.50	-9.25		-9.50		-9.50		-9.25		-9.50		-9.25		-9.50		-9.25		-9.50		-9.50	
-9.75	-9.50		-9.75		-9.75		-9.50		-9.75		-9.50		-9.75		-9.50		-9.75		-9.75	
-10.00	-9.75		-10.00		-10.00		-9.75		-10.00		-9.75		-10.00		-9.75		-10.00		-10.00	
-10.25	-10.00		-10.25		-10.25		-10.00		-10.25		-10.00		-10.25		-10.00		-10.25		-10.25	
-10.50	-10.25		-10.50		-10.50		-10.25		-10.50		-10.25		-10.50		-10.25		-10.50		-10.50	
-10.75	-10.50		-10.75		-10.75		-10.50		-10.75		-10.50		-10.75		-10.50		-10.75		-10.75	
-11.00	-10.75		-11.00		-11.00		-10.75		-11.00		-10.75		-11.00		-10.75		-11.00		-11.00	
-11.25	-11.00		-11.25		-11.25		-11.00		-11.25		-11.00		-11.25		-11.00		-11.25		-11.25	
-11.50	-11.25		-11.50		-11.50		-11.25		-11.50		-11.25		-11.50		-11.25		-11.50		-11.50	
-11.75	-11.50		-11.75		-11.75		-11.50		-11.75		-11.50		-11.75		-11.50		-11.75		-11.75	
-12.00	-11.75		-12.00		-12.00		-11.75		-12.00		-11.75		-12.00		-11.75		-12.00		-12.00	-1.75

※この表は各眼鏡度数におけるトーリックコンタクトレンズのトライアルの目安を示したものです(頂点間距離は12mmを想定)。トライアルの上、見え方を確認し、最終度数を決定してください。一部、製品の製作度数範囲に含まれていない部分があります。

MyDay[®] toric

マイデイ[®] トーリック



販売名:マイデイ 承認番号:22700BZX00320000

Biofinity[®] toric



バイオフィニティ[®] トーリック



販売名:バイオフィニティ 承認番号22200BZX00714A01

製作範囲表 (乱視用ソフトコンタクトレンズ)

販売用 / トライアル

		1day			2week		
製品	 マイデイ® トーリック			 バイオフィニティ® トーリック			
円柱度数	-0.75 -1.25 -1.75 -2.25	-0.75 -1.25 -1.75	-2.25	-0.75 -1.25 -1.75 -2.25	-0.75 -1.25 -1.75	-2.25	
円柱軸	90° 180°	10°20°70°80°90°100° 110°160°170°180°	10°20°90° 160°170°180°	90° 180°	10°20°90° 160°170°180°		
球面度数	+6.00						
	+5.75						
	+5.50						
	+5.25						
	+5.00						
	+4.75						
	+4.50						
	+4.25						
	+4.00						
	+3.75						
	+3.50						
	+3.25						
	+3.00						
	+2.75						
	+2.50						
	+2.25						
	+2.00						
	+1.75						
	+1.50						
	+1.25						
	+1.00						
	+0.75						
	+0.50						
	+0.25						
	±0.00						
	-0.25						
	-0.50						
	-0.75						
	-1.00						
	-1.25						
	-1.50						
	-1.75						
	-2.00						
	-2.25						
	-2.50						
	-2.75						
	-3.00						
	-3.25						
	-3.50						
	-3.75						
	-4.00						
	-4.25						
	-4.50						
	-4.75						
	-5.00						
	-5.25						
	-5.50						
	-5.75						
	-6.00						
	-6.50						
	-7.00						
	-7.50						
	-8.00						
	-8.50						
	-9.00						
	-9.50						
	-10.00						